




UBA52 Antibody

Product Code	CSB-PA035362
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62987
Immunogen	Synthesized peptide derived from internal of Human RPL40.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total RPL40 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	60S ribosomal protein L40; CEP52;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBA52

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using RPL40 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.