




PGD Antibody

Product Code	CSB-PA787749
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52209
Immunogen	Synthesized peptide derived from internal of Human PGD.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PGD protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	6-phosphogluconate dehydrogenase; 6-phosphogluconate dehydrogenase; decarboxylating; 6PGD; PGDH
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PGD

Image



Western blot analysis of extracts from A549 cells, using PGD antibody.