





PGLS Antibody

Product Code	CSB-PA222487
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95336
Immunogen	Synthesized peptide derived from internal of Human PGLS.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PGLS protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	6-phosphogluconolactonase; EC 3.1.1.31; 6PGL;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PGLS

Image



Western blot analysis of extracts from COLO cells, using PGLS antibody.



Immunohistochemistry analysis of paraffin-embedded human lymph node tissue, using PGLS antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.