




ABHD11 Antibody

Product Code	CSB-PA973412
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFV4
Immunogen	Synthesized peptide derived from internal of Human ABHD11.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABHD11 protein
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABHDB; abhydrolase domain containing 11; EC 3.-.-.-; PP1226; WBSCR21
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD11

Image



Western blot analysis of extracts from HUVEC cells and COLO cells, using ABHD11 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.