




ABHD12 Antibody

Product Code	CSB-PA589699
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2K0
Immunogen	Synthesized peptide derived from internal of Human ABHD12.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ABHD12 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABD12; ABHD12; abhydrolase domain containing 12;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD12

Image



Western blot analysis of extracts from COLO cells, using ABHD12 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.