




ABHD4 Antibody

Product Code	CSB-PA901885
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TB40
Immunogen	Synthesized peptide derived from internal of Human ABHD4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ABHD4 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Abhydrolase domain-containing protein 4; EC 3.1.1.-; Lyso-N-acylphosphatidylethanolamine lipase; Alpha/beta-hydrolase 4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD4

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ABHD4 antibody.