




ABHD6 Antibody

Product Code	CSB-PA449719
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV23
Immunogen	Synthesized peptide derived from internal of Human ABHD6.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ABHD6 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Abhydrolase domain-containing protein 6; EC 3.-.-.-; ABHD6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD6

Image



Western blot analysis of extracts from COLO cells, using ABHD6 antibody.



Immunofluorescence analysis of A549 cells, using ABHD6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.