




ACBD6 Antibody

Product Code	CSB-PA050508
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BR61
Immunogen	Synthesized peptide derived from internal of Human ACBD6.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ACBD6 protein
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	acyl-CoA-binding domain-containing 6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACBD6

Image



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using ACBD6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.