





ACOT1 Antibody

Product Code	CSB-PA132650
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TX2
Immunogen	Synthesized peptide derived from internal of Human ACOT1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ACOT1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ACH2; CTE-1; EC 3.1.2.2; LACH2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACOT1

Image



Western blot analysis of extracts from Jurkat cells, using ACOT1 antibody.



Immunofluorescence analysis of MCF-7 cells, using ACOT1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.