




# ACOT12 Antibody

<b>Product Code</b>	CSB-PA066461
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WYK0
<b>Immunogen</b>	Synthesized peptide derived from internal of Human ACOT12.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total ACOT12 protein.
<b>Tested Applications</b>	ELISA, IHC, IF; IHC: 1:50-1:100, IF: 1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Acyl-coenzyme A thioesterase 12; EC 3.1.2.1; Acyl-CoA thioesterase 12; Acyl-CoA thioester hydrolase 12; Cytoplasmic acetyl-CoA hydrolase 1
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ACOT12

## Image



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using ACOT12 antibody.



Immunofluorescence analysis of A549 cells, using ACOT12 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.