




ACOT12 Antibody

Product Code	CSB-PA066461
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYK0
Immunogen	Synthesized peptide derived from internal of Human ACOT12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ACOT12 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Acyl-coenzyme A thioesterase 12; EC 3.1.2.1; Acyl-CoA thioesterase 12; Acyl-CoA thioester hydrolase 12; Cytoplasmic acetyl-CoA hydrolase 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACOT12

Image



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using ACOT12 antibody.



Immunofluorescence analysis of A549 cells,
using ACOT12 antibody.