




ACOT2 Antibody

Product Code	CSB-PA570173
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49753
Immunogen	Synthesized peptide derived from internal of Human ACOT2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ACOT2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	acyl-coenzyme A thioesterase 2; CTE-Ia; EC 3.1.2.2; peroxisomal acyl-coenzyme A thioester hydrolase 2a; peroxisomal long-chain acyl-coA thioesterase 2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACOT2

Image



Western blot analysis of extracts from Jurkat cells and HuvEc cells, using ACOT2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.