





# ACOT4 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA051355   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q8N9L9   |
| <b>Immunogen</b>           | Synthesized peptide derived from internal of Human ACOT4.  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Specificity</b>         | The antibody detects endogenous levels of total ACOT4 protein.   |
| <b>Tested Applications</b> | ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500   |
| <b>Form</b>                | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                       |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | acyl-CoA thioesterase 4; acyl-coenzyme A thioesterase 4; EC 3.1.2.2; peroxisomal acyl coenzyme A thioester hydrolase Ib; peroxisomal long-chain acyl-CoA thioesterase Ib |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | ACOT4  |

## Image



Western blot analysis of extracts from HeLa cells and COLO cells, using ACOT4 antibody.



Immunofluorescence analysis of MCF-7 cells,  
using ACOT4 antibody.