





ACOT4 Antibody

Product Code	CSB-PA051355
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9L9
Immunogen	Synthesized peptide derived from internal of Human ACOT4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ACOT4 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	acyl-CoA thioesterase 4; acyl-coenzyme A thioesterase 4; EC 3.1.2.2; peroxisomal acyl coenzyme A thioester hydrolase Ib; peroxisomal long-chain acyl-CoA thioesterase Ib
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACOT4

Image



Western blot analysis of extracts from HeLa cells and COLO cells, using ACOT4 antibody.



Immunofluorescence analysis of MCF-7 cells, using ACOT4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.