




ADAM32 Antibody

Product Code	CSB-PA552451
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C14296
Immunogen	Synthesized peptide derived from Internal of Human ADAM32.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ADAM32 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	a disintegrin and metalloprotease domain 32; a disintegrin and metalloproteinase domain 32; ADA32; ADAM 32; ADAM metalloproteinase domain 32
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADAM32

Image



Western blot analysis of extracts from 293 cells, using ADAM32 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.