




ADAMTS18 Antibody

Product Code	CSB-PA999160
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE60
Immunogen	Synthesized peptide derived from internal of Human ADAMTS18.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ADAMTS18 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	A disintegrin and metalloproteinase with thrombospondin motifs 18; ADAM-TS 18; ADAM-TS18; ADAMTS-18EC=3.4.24.-;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADAMTS18

Image



Western blot analysis of extracts from 3T3 cells, using ADAMTS18 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.