




BST1 Antibody

Product Code	CSB-PA052030
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q10588
Immunogen	Synthesized peptide derived from internal of Human BST1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total BST1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADP-ribosyl cyclase 2; bone marrow stromal antigen 1; BST-1; cADPr hydrolase 2; CD157 antigen
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BST1

Image



Western blot analysis of extracts from Jurkat cells and COS cells, using BST1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.