




AKR1A1 Antibody

Product Code	CSB-PA967994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14550
Immunogen	Synthesized peptide derived from C-terminal of Human AKR1A1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AKR1A1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AK1A1; AKR1A1; alcohol dehydrogenase; Alcohol dehydrogenase [NADP+]; Aldehyde reductase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKR1A1

Image



Western blot analysis of extracts from COLO205 cells and NIH-3T3 cells, using AKR1A1 antibody.