




ADH7 Antibody

Product Code	CSB-PA098332
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40394
Immunogen	Synthesized peptide derived from internal of Human ADH7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ADH7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADH-4; alcohol dehydrogenase class 4 mu/sigma; alcohol dehydrogenase class IV mu/sigma; EC 1.1.1.1; gastric alcohol dehydrogenase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADH7

Image



Western blot analysis of extracts from COS cells, using ADH7 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.