




AKR1CL1 Antibody

Product Code	CSB-PA175997
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T2L2
Immunogen	Synthesized peptide derived from internal of Human AKR1CL1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total AKR1CL1 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Aldo-keto reductase family 1 member C-like protein 1; EC 1.1.1.-;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKR1C8P

Image



Immunofluorescence analysis of A549 cells, using AKR1CL1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.