




ACER2 Antibody

Product Code	CSB-PA293630
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5QJU3
Immunogen	Synthesized peptide derived from internal of Human ASAH3L.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total ASAH3L protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Alkaline ceramidase 2; AlkCDase 2; Alkaline CDase 2; haCER2; EC=3.5.1.23
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACER2

Image



Western blot analysis of extracts from COS7 cells, using ASAH3L antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.