




SGCA Antibody

Product Code	CSB-PA993648
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16586
Immunogen	Synthesized peptide derived from internal of Human SGCA.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SGCA protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Alpha-sarcoglycan; Alpha-SG; Adhalin; 50 kDa dystrophin-associated glycoprotein; 50DAG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SGCA

Image



Western blot analysis of extracts from Jurkat cells, HepG2 cells and HeLa cells, using SGCA antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.