




# NECAB3 Antibody

<b>Product Code</b>	CSB-PA148617
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96P71
<b>Immunogen</b>	Synthesized peptide derived from internal of Human NECAB3.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total NECAB3 protein.
<b>Tested Applications</b>	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	AB2BP; amyloid beta A4 protein-binding family A member 2-binding; APBA2BP; dJ63M2.4; dJ63M2.5
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NECAB3

## Image



Western blot analysis of extracts from COLO cells, HepG2 cells and HT-29 cells, using NECAB3 antibody.



Immunofluorescence analysis of A549 cells, using NECAB3 antibody.

### Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.