




ANGPTL7 Antibody

Product Code	CSB-PA028244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43827
Immunogen	Synthesized peptide derived from C-terminal of Human ANGPTL7.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ANGPTL7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	angiopoietin-like 7; ANGL7; AngX; CDT6; cornea-derived 6
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANGPTL7

Image



Western blot analysis of extracts from HepG2 cells, Jurkat cells, COLO cells and HUVEC cells, using ANGPTL7 antibody.