




# ANGPTL7 Antibody

|                            |                                                                                                                                                    |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA028244                                                                                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                      |
| <b>Uniprot No.</b>         | O43827                                                                                                                                             |
| <b>Immunogen</b>           | Synthesized peptide derived from C-terminal of Human ANGPTL7.                                                                                      |
| <b>Raised In</b>           | Rabbit                                                                                                                                             |
| <b>Species Reactivity</b>  | Human, Mouse                                                                                                                                       |
| <b>Specificity</b>         | The antibody detects endogenous levels of total ANGPTL7 protein.                                                                                   |
| <b>Tested Applications</b> | ELISA, WB; WB: 1:500-1:3000                                                                                                                        |
| <b>Form</b>                | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| <b>Clonality</b>           | Polyclonal                                                                                                                                         |
| <b>Alias</b>               | angiopoietin-like 7; ANGL7; AngX; CDT6; cornea-derived 6                                                                                           |
| <b>Product Type</b>        | Polyclonal Antibody                                                                                                                                |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                               |
| <b>Target Names</b>        | ANGPTL7                                                                                                                                            |

## Image



Western blot analysis of extracts from HepG2 cells, Jurkat cells, COLO cells and HUVEC cells, using ANGPTL7 antibody.

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.