





ANKRD20A3 Antibody

Product Code	CSB-PA136100
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VUR7
Immunogen	Synthesized peptide derived from internal of Human ANKRD20A3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ANKRD20A3 protein
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ankyrin repeat domain-containing protein 20A3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD20A3

Image



Western blot analysis of extracts from 3T3 cells, using ANKRD20A3 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ANKRD20A3 antibody.