





KANK2 Antibody

Product Code	CSB-PA178754
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63ZY3
Immunogen	Synthesized peptide derived from internal of Human KANK2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total KANK2 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ankyrin repeat domain-containing protein 25; SRC-1-interacting protein; SRC1-interacting protein; Matrix-remodeling-associated protein 3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KANK2

Image



Immunohistochemistry analysis of paraffin-embedded human heart tissue using KANK2 antibody.



Immunofluorescence analysis of A549 cells, using KANK2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.