




APOA5 Antibody

Product Code	CSB-PA206029
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Q788
Immunogen	Synthesized peptide derived from internal of Human APOA5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total APOA5 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Apolipoprotein A-V; Apo-AV; ApoA-V; Apolipoprotein A5; Regeneration-associated protein 3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APOA5

Image



Western blot analysis of extracts from 293 cells, using APOA5 antibody.