




APOL1 Antibody

Product Code	CSB-PA923087
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14791
Immunogen	Synthesized peptide derived from internal of Human APOL1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total APOL1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Apo-L; APOL; ApoL-I; Apolipoprotein L; Apolipoprotein L-I
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APOL1

Image



Western blot analysis of extracts from A549 cells, using APOL1 antibody.