




AQP1 Antibody

Product Code	CSB-PA118586
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29972
Immunogen	Synthesized peptide derived from internal of Human AQP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total AQP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AQP-1; AQP-CHIP; aquaporin 1; aquaporin-1; aquaporin-CHIP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AQP1

Image



Western blot analysis of extracts from rat brain cells, using AQP1 antibody.