





ARSD Antibody

Product Code	CSB-PA215417
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51689
Immunogen	Synthesized peptide derived from internal of Human ARSD.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ARSD protein
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Arylsulfatase D; EC 3.1.6.-; ASD;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARSD

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using ARSD antibody.



Immunofluorescence analysis of MCF-7 cells, using ARSD antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.