




ATP5H Antibody

Product Code	CSB-PA146706
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75947
Immunogen	Synthesized peptide derived from internal of HumanATP5H.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ATP5H protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP synthase D chain; mitochondrial; ATP synthase; H ⁺ transporting; mitochondrial F ₀ complex
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP5H

Image



Western blot analysis of extracts from HepG2 cells and Jurkat cells, using ATP5H antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.