




ATP5I Antibody

Product Code	CSB-PA796067
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56385
Immunogen	Synthesized peptide derived from C-terminal of Human ATP5I.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ATP5I protein
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP synthase e chain; mitochondrial; ATP5K; EC 3.6.3.14; F1F0-ATP synthase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP5I

Image



Western blot analysis of extracts from COLO cells, using ATP5I antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.