




ABCA6 Antibody

Product Code	CSB-PA988363
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N139
Immunogen	Synthesized peptide derived from C-terminal of Human ABCA6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABCA6 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABC transporter ABCA6; ABCA6; ATP-binding cassette A6; ATP-binding cassette sub-family A member 6; ATP-binding cassette
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCA6

Image



Western blot analysis of extracts from 293 cells, using ABCA6 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.