




ABCA8 Antibody

Product Code	CSB-PA847031
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94911
Immunogen	Synthesized peptide derived from internal of Human ABCA8.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ABCA8 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family A member 8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCA8

Image



Western blot analysis of extracts from HUVEC cells, using ABCA8 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.