





ABCB7 Antibody

Product Code	CSB-PA575138
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75027
Immunogen	Synthesized peptide derived from internal of Human ABCB7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABCB7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABC transporter 7 protein; ABC7; Anemia; sideroblastic; with spinocerebellar ataxia
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCB7

Image



Western blot analysis of extracts from HepG2 cells, using ABCB7 antibody.



Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using ABCB7 antibody.