




ABCD4 Antibody

Product Code	CSB-PA970746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14678
Immunogen	Synthesized peptide derived from internal of Human ABCD4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABCD4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family D member 4; Peroxisomal membrane protein 69; PMP69; Peroxisomal membrane protein 1-like; PXMP1-L
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCD4

Image



Western blot analysis of extracts from HepG2 cells, using ABCD4 antibody.