




ABCF1 Antibody

Product Code	CSB-PA554269
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE71
Immunogen	Synthesized peptide derived from C-terminal of Human ABCF1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous levels of total ABCF1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ABC27; ABC50; ABCF1; ATP-binding cassette 50; ATP-binding cassette 50 (TNF-alpha stimulated)
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCF1

Image



Western blot analysis of extracts from HeLa/COLO205 cells, using ABCF1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.