




DDX3Y Antibody

Product Code	CSB-PA947616
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15523
Immunogen	Synthesized peptide derived from internal of Human DDX3Y.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total DDX3Y protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase DDX3Y; DBY; DDX; DEAD (Aspartate-glutamate-alanine-aspartate) box polypeptide; Y chromosome
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DDX3Y

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using DDX3Y antibody.