




DDX51 Antibody

Product Code	CSB-PA053041
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8A6
Immunogen	Synthesized peptide derived from C-terminal of Human DDX51.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DDX51 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase DDX51; DEAD box protein 51; EC 3.6.1.-;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DDX51

Image



Western blot analysis of extracts from HepG2 cells, using DDX51 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.