




BAGE3 Antibody

Product Code	CSB-PA588687
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y29
Immunogen	Synthesized peptide derived from internal of Human BAGE3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BAGE3 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cancer/testis antigen 2.3; CT2.3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAGE3

Image



Immunohistochemistry analysis of paraffin-embedded human testis tissue using BAGE3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.