




GCNT3 Antibody

Product Code	CSB-PA976989
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95395
Immunogen	Synthesized peptide derived from internal of Human GCNT3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GCNT3 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1; 6-N-acetylglucosaminyltransferase 3; EC 2.4.1.102; EC 2.4.1.150
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GCNT3

Image



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using GCNT3 antibody.



Immunofluorescence analysis of A549 cells, using GCNT3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.