





GCNT7 Antibody

Product Code	CSB-PA596794
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZNI0
Immunogen	Synthesized peptide derived from internal of Human GCNT7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GCNT7 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1; 6-N-acetylglucosaminyltransferase 7; EC 2.4.1.-; C20orf105
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GCNT7


Image



Western blot analysis of extracts from K562 cells, using GCNT7 antibody.



Immunohistochemistry analysis of paraffin-embedded human ovary tissue using GCNT7 antibody.



Immunofluorescence analysis of MCF-7 cells, using GCNT7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.