





B3GALTL Antibody

Product Code	CSB-PA786105
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Y288
Immunogen	Synthesized peptide derived from C-terminal of Human B3GALTL.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total B3GALTL protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	B3Glc-T; B3GLT; B3GTL; beta 1; 3-galactosyltransferase-like
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	B3GLCT

Image



Western blot analysis of extracts from HepG2 cells, Jurkat cells, MCF-7 cells and HUVEC cells, using B3GALTL antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using B3GALTL antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.