




GLB1 Antibody

Product Code	CSB-PA146106
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16278
Immunogen	Synthesized peptide derived from internal of Human GLB1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GLB1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Beta-galactosidase; EC=3.2.1.23; Acid beta-galactosidase; Lactase; Elastin receptor 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GLB1

Image



Western blot analysis of extracts from 3T3 cells, using GLB1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.