




BMP10 Antibody

Product Code	CSB-PA198870
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95393
Immunogen	Synthesized peptide derived from Internal of Human BMP10.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BMP10 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BMP-10; BMP10; Bone morphogenetic protein 10; MGC126783;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BMP10

Image



Western blot analysis of extracts from Jurkat cells, using BMP10 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.