





HGH1 Antibody

Product Code	CSB-PA989731
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTY7
Immunogen	Synthesized peptide derived from internal of Human BRP16.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BRP16 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Brain protein 16;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HGH1

Image



Western blot analysis of extracts from COLO cells, using BRP16 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using BRP16 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.