




MPC2 Antibody

Product Code	CSB-PA067696
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95563
Immunogen	Synthesized peptide derived from internal of Human BRP44.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total BRP44 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Brain protein 44;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MPC2

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using BRP44 antibody.