




BAIAP2L1 Antibody

Product Code	CSB-PA794497
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHR4
Immunogen	Synthesized peptide derived from internal of Human BAIAP2L1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BAIAP2L1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAIAP2L1; BI2L1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAIAP2L1

Image



Western blot analysis of extracts from HepG2 cells, using BAIAP2L1 antibody.