




BAIAP2L2 Antibody

Product Code	CSB-PA000226
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXY1
Immunogen	Synthesized peptide derived from internal of Human BAIAP2L2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total BAIAP2L2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAI1-associated protein 2-like 2; BI2L2; brain-specific angiogenesis inhibitor 1-associated 2-like 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAIAP2L2

Image



Western blot analysis of extracts from RAW264.7 cells, using BAIAP2L2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.